Multifunction System **Kit**

USER MANUAL

MANUAL DE USO / GUIDE DE L'UTILISATEUR / MANUALE D'USO /
MANUAL DE UTILIZAÇÃO / GEBRAUCHSANWEISUNG FÜR /
РУКОВОДСТВО ПО ИСПОЛЬЗОВАНИЮ НАБОРА

| | | | | | | | | |
|------------------------------|---------|------|------|-------|-------|-------|-------|-------|
| NS Nose Stimulator | MM SIZE | 00 | 0 | 1 | 2 | 3 | | |
| | | 8,81 | 9,51 | 11,02 | 12,02 | 13,03 | 14,03 | 15,03 |

Estimulador nasal / Dilatateur nasal / Dispositivo nasale / Estimulador nasal / Nasen stimulator / Носовой стимулятор

| | | | | | | | | |
|----------------------------------|---------|--|------|-------|-------|-------|--------|--------|
| PO Permeable Obturator | MM SIZE | | 1 | 2 | 3 | 4 | | |
| | | | 68,3 | 77,22 | 86,79 | 94,68 | 107,04 | 118,79 |

Obturador permeable / Obturateur perméable / Otturatore permeabile / Obturador permeável / Durchlässige verschlussplatte / Проницаемый obturator

| | | | | | | | | |
|-----------------------------|---------|--|------|-------|-------|-------|--------|--------|
| LS Lip Stimulator | MM SIZE | | 1 | 2 | 3 | 4 | | |
| | | | 68,3 | 77,22 | 86,79 | 94,68 | 107,04 | 118,79 |

Estimulador labial / Stimulateur labial / Stimolatore labiale / Estimulador labial / Lippenstimulator / Губной стимулятор

| | | | | | | | | |
|------------------------------|---------|--|----|------|-------|-------|-------|-----|
| MR Muscle Relaxant | MM SIZE | | 1 | 2 | 3 | 4 | | |
| | | | 83 | 93,6 | 104,2 | 114,8 | 125,4 | 136 |

Relajante muscular / Relaxant musculaire / Dispositivo per rilassare i muscoli / Relaxante muscular / Muskelentspanner / Приспособление для расслабления мышц

| | | | | | | | | |
|-------------------------|---------|--|------|-------|-------|-------|-------|----|
| BP Bite Plate | MM SIZE | | 1 | 2 | 3 | 4 | | |
| | | | 28,8 | 32,64 | 36,48 | 40,32 | 44,16 | 48 |

Plano de mordida / Plan de morsure / Morso / Placa de mordida / Bissfläche / Накусочная площадка

| | | | | | | | | |
|------------------------|---------|--|-------|-------|------|-------|-------|-------|
| OB Open Bite | MM SIZE | | 1 | 2 | 3 | 4 | | |
| | | | 55,24 | 58,92 | 62,6 | 66,29 | 69,97 | 73,65 |

Mordida abierta / Ouvert morsure / Bite aperto / Mordida aberta / Offene bisschiene / Открытый прикус




INSTRUCTIONS FOR USE

INSTRUCCIONES DE USO / COMMENT L'UTILISER / COME SI USA /
MÉTODO DE UTILIZAÇÃO / ANWENDUNG / СПОСОБ ПРИМЕНЕНИЯ

NOSE STIMULATOR

Estimulador nasal / Dilatateur nasal / Dispositivo nasale / Estimulador nasal / Nasen stimulator / Носовой стимулятор

Fig. 1

ORAL OBTURATORS



Obturador / Obturateur / Otturatore / Obturador / Verschlussplatte / Обтurator

Fig. 2

1. Permeable / Permeable /
Perméable / Permeabile / Permeável /
Durchlässige / Проницаемый




Fig. 3

LIP STIMULATOR

Estimulador labial / Stimulateur labial / Stimolatore labiale / Estimulador labial / Lippenstimulator / Губной стимулятор


Fig. 4

MUSCLE RELAXANT

Relajante muscular / Relaxant musculaire / Dispositivo per rilassare i muscoli / Relaxante muscular / Muskelentspanner / Приспособление для расслабления мышц


Fig. 5

BITE PLATE

Plano de mordida / Plan de morsure / Morso / Placa de mordida / Bissfläche / Накусочная площадка


Fig. 6

ANTI-BRUXISM DEVICE

Dispositivo Antibruxista / Dispositif Anti-bruxisme / Dispositivo Anti-bruxismo / Dispositivo Antibruxismo / Bruxomanievorrichtung / Аппарат Против Бруксизма


Fig. 7

OPEN BITE

Mordida abierta / Ouvert morsure / Bite aperto / Mordida aberta / Offene bissschiene / Открытый прикус

Fig. 8

Multifunction System Kit

ENGLISH

INTRODUCTION

The “programmed stimulotherapy” protocol, resulting from the evaluation of oral function, based on the objective measurement and quantification of collapsed nares, the degree of adenoid and tonsillar hypertrophy, lingual mobility, and patterns of deglutition, involves the use of any of the MFS devices in a prescribed sequence of automated exercises with the stated objective of promoting the normalisation of the triumvirate of oral function: respiration, deglutition, and masticatory pattern.

Clinical observation of the relevance of oral motor functions led to the idea of creating muscular stimuli as a basis for the treatment of impaired oral functions.

These devices serve as a basis for an oral motor performance rehabilitation programme, both in patients with removable or fixed oral devices, or in patients not wearing any type of oral device. This is particularly helpful in pediatric patients, since it contributes to the total oral health programme. In addition, these devices may be clinically used as a coadjuvant for specific anti-snoring, sleep apnea and bruxism treatments.

The MFS kit offers a wide range of therapeutic tools to the orthodontic clinical practice, providing a new clinical approach and new treatment perspectives. The kit may be used in early diagnostics and treatment and even in late stage treatment. It also constitutes a new method for oral disease prevention in clinical orthodontics.

The MFS KIT features 8 distinct devices, seven of which come in 6 different sizes and one in 7 sizes. The kit includes 2 units per size.

All MULTIFUNCTION SYSTEM Kit devices are class 1 sanitary products, in compliance with the EU REGULATIONS governing sanitary products.

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NOSE STIMULATORS

1. WHAT IS A NOSE STIMULATOR?

The nose stimulator consists of two tubes that are joined by a stabilizing tape, and comprises a flat surface that makes direct contact with the nose septum, an external convex surface portion that tenses the nasal alae, a tab that stimulates the muscular insertions in the nasal alar area, and an external rim to prevent the tubes from accidentally hitting against the nose. The device comes in seven sizes, two for children (00 and 01) and five for adolescents and adults (1, 2, 3, 4, and 5).

It stimulates perinasal muscular insertions in the alar region of the nose, reduces nasal airflow resistance, enhancing nasal airflow, and restructures the nasal cartilage, thereby harmonizing the shape of the nasal pyramid. It is used to promote nasal airflow in patients who present with:

- Nasal passages collapse during inhalation
- Mouth breathing
- Snoring problems
- Sleep apnea

2. INDICATIONS

The nose stimulator is especially indicated for patients with nasal alar collapse resulting in restriction of airflow through the nasal orifice. It is moreover indicated for mouth-breathing subjects to promote nasal breathing. Given that the device improves nasal airflow, it is also indicated for snoring subjects to palliate the symptoms, since enhanced nasal airflow increases the oxygen available to the lung alveoli.

Only a licensed dentist may perform a complete evaluation of each case history and prescribe the use of this device, whether it is used alone or in combination with other devices.

3. WHEN AND HOW TO USE A NOSE STIMULATOR

Guidelines to ensure that the nose stimulator is correctly used:

1. Select the adequate nose stimulator size, measuring the diameter of the dilated nostrils using an orthodontic Vernier caliper.
2. The diameter of the nose stimulator must be equal to the diameter of the nostrils. In case of doubt, select a device of smaller diameter.
3. Coat the nasal cavities with petroleum jelly before introducing the chosen nose stimulator to facilitate introduction.
4. A night-time regimen is indicated.
5. After several months of use, the patient is given a larger size. The size of the device will be progressively increased throughout the treatment until maximum nostril dilatation is obtained.

Instructions for Use: refer to page 0, fig. 1

4. RESULTS

1. The immediate result is increased nasal airflow volume, which in turn results in an improvement in:
 - a. Nasal breathing
 - b. Snoring and sleep apnea
2. The long-term result (7-9 months) is an improvement in nostril dilatation (diameter) during inhalation, thanks to the improved performance of the perinasal muscle group after stimulation of motor units with the device.
3. After nine months, the effects of the nose stimulators may wane, in which case the following regimen will be indicated:
 - a. Snoring subjects, especially those patients who also suffer from sleep apnea, should continue treatment as advised.
 - b. Mouth-breathing subjects should use the device every six months –for a three-month period- until the breathing problem is fully normalised.
4. Nose stimulators may provide a partial or total solution to sleep apnea patients. In some cases, it is convenient to use the device in combination with a mandibular propulsion appliance to advance the tongue.

5. SIZES

Nose stimulators come in five sizes for adult subjects and two for pediatric patients over 5 years of age. The devices are numbered according to size: Sizes 00 and 0 –for pediatric patients- and sizes 1, 2, 3, 4, and 5 for patients over 12 years of age.

Recommended size: Initially, each patient will be given a size smaller than nostril diameter. The size of the device will be progressively increased throughout the treatment in order to ensure a steady adaptation of the stimulator to nostril diameter while obtaining gradual, more comfortable effects on the patient's nose.

This device may be used in combination with fixed orthodontic appliances.

6. CONTRAINDICATIONS

Should not be used in patients under five years of age.

7. WARNINGS

Nose stimulators should not be placed in direct contact with any heat source, otherwise the devices could be deformed and their elasticity impaired.

8. OBSERVATIONS

The results confirm that nose stimulators, designed for night-time use, have clinically important effects. They are also an effective therapeutic tool for functional respiratory re-education to correct mouth breathing and for nasal alar collapse correction, thanks to its immediate effects and to its medium-term action.

ORAL OBTURATORS

1. WHAT IS AN ORAL OBTURATOR?

The oral obturator is used to obturate the oral cavity and gradually impede mouth breathing, which may cause malocclusion as well as relapse after orthodontic treatment. The device is a simple shield that progressively impedes the flow of air through the oral cavity, thereby aiding mouth-breathing patients.

It has been especially designed to adapt itself to the shape of the dental arches and works as a seal that prevents mouth breathing. The oral obturator features a peripheral rim portion –extending around the upper and lower end- to allow the patient to exercise the labial muscles, which are very important for normalisation of nasal breathing.

There are three types of obturators -permeable, semi-permeable and impermeable- depending on whether they include large breathing holes, small holes, or no holes at all. Each type of obturator comes in six sizes.

- Permeable oral obturator, featuring large breathing holes to allow limited air flow.
- Semi-permeable oral obturator, including smaller breathing holes to control and limit and limit air flow through the mouth even further.
- Impermeable oral obturator, with no breathing holes to totally impede air flow through the mouth.

2. INDICATIONS

Oral obturators are indicated for mouth-breathing patients, to aid in normalising the breathing problem by means of progressive respiratory re-education to correct mouth breathing and promote nasal route breathing.

Only a licensed dentist may perform a complete evaluation of each case history and prescribe the use of this device, whether it is used alone or in combination with other devices.

3. WHEN AND HOW TO USE AN ORAL OBTURATOR

Oral obturators are used in combination with nasal stimulators to avoid producing a “suffocating” effect in the patient. Their use is also recommended once the physiological causes of mouth breathing have been eliminated: rhinitis, nasal septum deviation, hypertrophic adenoids or tonsils, among the most frequent causes.

Guidelines to ensure that the oral obturator is correctly used:

1. The dental arch perimeter between the first permanent molars measurement is taken from the patient’s mouth or from the mold of the patient’s dental configuration.
2. The adequate oral obturator size is chosen based on this measurement. In case of doubt, select a smaller device than your measured size.

3. Use a permeable obturator to initiate the respiratory pattern normalisation process. It is essential to combine the use of the obturator with the nasal stimulator.
4. The obturator is placed by the patient in the buccal vestibule, between the dental arches, and the lips and cheeks. A night-time regimen is indicated, and if the patient experiences discomfort when using the device at first, the patient should be instructed to wear it at home during the day to get used to it.
5. If the patient has already adapted to using the permeable obturator, after three or four months the patient is instructed to use the semi-permeable obturator.
6. After three or four months of using the semi-permeable obturator, the patient is ready to use the impermeable device.
7. The improvement in labial competence will be also clinically evaluated, and if full recovery of labial competence is not achieved, then the lip stimulator must be used.
8. This medical device may be used in combination with fixed orthodontic appliances.

Types of obturators: refer to page 0, fig. 2

Instructions for Use: refer to page 0, fig. 3

4. SIZES

Oral obturators come in 6 sizes (numbered from 1 to 6): the three first sizes are typically used in the pediatric population, while the remaining sizes have been designed for use in adult patients.

Since oral obturators work passively, that is to say, without exerting any force since its main objective is to gradually reduce the passage of air through the mouth, it is very important to select the right size.

5. CONTRAINDICATIONS

Should not be used in patients under five years of age, or in mouth-breathing patients who still suffer from nasal blockage, leading to breathing problems.

6. WARNINGS

Oral obturators should not be placed in direct contact with any heat source, otherwise the devices could be deformed and their elasticity impaired.

7. OBSERVATIONS

Oral obturators has a double function: to control the passage of air through the mouth and to stimulate the labial muscles. These effects are gradually achieved and a reasonable amount of time must elapse before the results may be clinically confirmed.

LIP STIMULATOR

1. WHAT IS A LIP STIMULATOR?

The lip stimulator is an oval-shaped device comprising an upper and lower front band. It acts in the back buccal vestibule area stimulating the orbicular muscles of the mouth. This stimulation causes the lips to draw together, which in turn produces the lengthening of the upper lip, thereby helping to correct labial incompetence. Indirectly, lip closing exercises cause the tongue to be posteriorly repositioned, an optimal tongue position to initiate the re-education of atypical deglutition.

2. INDICATIONS

Lip stimulators are indicated for patients with labial incompetence (parted lips at rest), with short upper lip (gingival smile), and/or with atypical deglutition involving forward tongue thrust and tongue pressure on the incisors.

3. WHEN AND HOW TO USE A LIP STIMULATOR

Guidelines to ensure that the lip stimulator is correctly used:

1. The dental arch perimeter between the first permanent molars measurement is taken from the patient's mouth or from the mold of the patient's dental configuration.
2. The adequate lip stimulator size is chosen based on this measurement. In case of doubt, select a smaller device than your measured size.
3. Once the lip stimulator is placed on the mouth, the patient will start exercising the lip muscles without being instructed by the health professional. This indicates that lip stimulation is working and the patient is able to achieve a competent lip seal.
4. A night-time regimen is indicated. Initially, and to get used to wearing the device, the patient may use it only during the day and at home.
5. To re-educate a low, forward tongue position the following points must be taken into account:
 - a. Remove any possible obstacle that may prevent tongue repositioning: lingual frenulum.
 - b. The lip stimulator is used until the tongue has repositioned itself backwards.
 - c. Bond several lingual buttons on to the surface of incisors or canines to stimulate the anterior palatal surface, and thereby promote tongue tip elevation. In some cases, a dental floss may be wound around

a button and extended towards a second button to obtain a more focused stimulation.

- d. After three months, bond several lingual buttons on to the surface of premolar teeth (stimulation of the posterior palatal surface) to help the dorsal surface of the tongue rise towards the palate during deglutition.

Instructions for Use: refer to page 0, fig. 4

4. RESULTS

1. The following effects resulting from the use of this device will be noticeable after four to six months of continuous use:
 - a. Competent lip seal (labial competence)
 - b. Lengthening of the upper lip
 - c. Backwards repositioning of the tongue. This is achieved through the lip closing exercises, performed with half-open mouth, which stimulate the muscles that reposition the tongue at rest to a normal position.

5. SIZES

The lip stimulator comes in six sizes, numbered from one to six. The smaller sizes are used in the pediatric population, while the larger sizes have been designed for use in adult patients.

The device's oval-shaped upper and lower front bands have a stimulating effect on the orbicular muscles due to their shape rather than to any force they exert. This is the reason why it is important to ensure that the correct stimulator size is being used.

6. CONTRAINDICATIONS

The use of lip stimulators in patients under five years of age is not recommended.

7. WARNINGS

Lip stimulators should not be placed in direct contact with any heat source, otherwise the devices could be deformed and their elasticity impaired.

8. OBSERVATIONS

Continuous use of the lip stimulator has reported relevant clinical results related to lip shape and position enhancement. The device is thereby highly indicated as a facial aesthetic treatment to increase lip muscle tone.

This medical device may be used in combination with fixed orthodontic appliances.

MUSCLE RELAXANT

1. WHAT IS A MUSCLE RELAXANT?

The muscle relaxant is a medical device included in the MFS line of health aid products. The device may be used alone or in combination with the bite plane.

The muscle relaxant has been especially designed to perform various functions in the mouth:

- It relaxes the masticatory muscles through the action of its lateral discs, which are vertically elevated to keep the mouth half open and the masseter muscle stretched and with less tone.
- It relaxes the buccinator muscle by pushing the cheeks outwards.
- It relaxes the orbicular muscles of the mouth by pushing the lips forward.

2. INDICATIONS

The muscle relaxant is especially indicated for patients with intense muscle contraction (brachycephalic) and in subject with peribuccal muscle contractions. Only a licensed dentist may perform a complete evaluation of each case history and prescribe the use of this device, whether it is used alone or in combination with other devices.

- Patients with intense muscle contraction (brachycephalic)
- Patients with tight lip muscles
- Patients with extremely tense perioral muscles

3. WHEN AND HOW TO USE A MUSCLE RELAXANT

Guidelines to ensure that the muscle relaxant is correctly used:

1. The dental arch perimeter between the first permanent molars measurement is taken from the patient's mouth or from the mold of the patient's dental configuration.
2. The adequate muscle relaxant size is chosen based on this measurement. In case of doubt, select a smaller device than your measured size.

3. This medical device must be placed in the patient's buccal vestibule, outside the dental arches.
4. A night-time regimen is indicated, although the patient may start wearing the muscle relaxant during daytime, to get used to the device.

Instructions for Use: refer to page 0, fig. 5

4. RESULTS

Clinical results are noticeable after six months of continuous use.

5. SIZES

The muscle relaxant comes in six sizes, numbered from one to six. The smaller sizes are used in the pediatric population, while the larger sizes have been designed for use in adult patients.

The dental arch perimeter between the first permanent superior molars measurement is taken from the patient's mouth and the adequate device size is chosen based on this measurement.

The device acts by applying a gentle pressure to the cheek area and along the lips to relax perioral muscles. This is the reason why it is important to ensure that the correct muscle relaxant size is being used.

6. CONTRAINDICATIONS

Should not be used in dolichocephalic patients (long narrow head) or in patients under five years of age.

7. WARNINGS

Muscle relaxants should not be placed in direct contact with any heat source, otherwise the devices could lose their physical characteristics.

8. OBSERVATIONS

This device may be used alone or in combination with the bite plane. The first option is useful for patients with moderately intense muscle contraction, whereas the second option is more suitable for patients presenting with extreme muscular contraction patterns and for patients with bruxism.

This device may be used in combination with fixed orthodontic appliances.

BITE PLATE

1. WHAT IS A BITE PLATE?

The bite plate is used only in combination with the muscle relaxant and is an integral part of the MFS anti-bruxism device. The device works as an anterior bite plate, an appliance regularly used in orthodonty, to raise the bite and to reduce deep overbite malocclusion of anterior teeth by extrusion of posterior teeth and posterior rotation of the jaw.

2. INDICATIONS

The bite plate is used only in combination with the muscle relaxant and is indicated in the palliative treatment of bruxism and in patients with deep anterior overbite.

Only a licensed dentist may perform a complete evaluation of each case history and prescribe the use of this device, whether it is used alone or in combination with other devices.

3. WHEN AND HOW TO USE A BITE PLATE

In combination with the muscle relaxant, the bite plate is an integral part of the MFS anti-bruxism device, and a night-time regimen is indicated. It is placed within the mouth in such a way that the muscle relaxant is located between the teeth, cheeks and lips, while the bite plate is inserted between the upper and lower incisors.

Instructions for Use: refer to page 0, fig. 6

4. SIZES

The bite plate comes in six sizes, numbered from one to six. The smaller sizes are used in the pediatric population, while the larger sizes have been designed for use in adult patients. Bite plate size matches the size of the muscle relaxant.

The bite plate works passively and is inserted between the upper and lower incisors.

To determine which size to use it is best to take the measurement from the mold of the patient's inferior dental arch.

5. CONTRAINDICATIONS

Should not be used in dolichocephalic patients (long narrow head) or in patients under five years of age.

6. WARNINGS

Bite plates should not be placed in direct contact with any heat source, otherwise the devices could be deformed and lose their physical characteristics.

7. OBSERVATIONS

The bite plate is used in combination with the muscle relaxant as a palliative treatment of bruxism. However, you should always consult a licensed dentist for a complete evaluation and diagnosis of the problem.

ANTI-BRUXISM DEVICE

1. WHAT IS THE ANTI-BRUXISM DEVICE?

The MFS anti-bruxism device combines the use of the muscle relaxant with the bite plate. These two devices are fixed together by the bite plate's mounting brackets which allow the device to be fixed to the anterior rectangular groove located at the muscle relaxant's anterior middle tape.

The following clinical results are obtained with the use of the MFS anti-bruxism device:

- Muscle relaxant-related results:
 - Relaxes the masticatory muscles through the action of its lateral discs, which are vertically elevated to keep the mouth half open and the masseter muscle stretched and with less tone.
 - It relaxes the buccinator muscle by pushing the cheeks outwards.
 - It relaxes the orbicular muscles of the mouth by pushing the lips forward.
- Bite plate-related results:
 - Reduces depth of anterior overbite
 - Generates extrusion of posterior teeth from occlusion
 - Creates posterior rotation of the jaw

2. INDICATIONS

The MFS anti-bruxism device is indicated in:

- Patients with bruxism
- Patients who have strong masticatory muscle activity patterns (brachycephalic)
- Patients with deep dental overbite. The device may be used in combination with fixed orthodontic appliances, and is much more comfortable to use than standard orthodontic appliances.

3. WHEN AND HOW TO USE THE ANTI-BRUXISM DEVICE

Guidelines to ensure that the MFS anti-bruxism device is correctly used:

1. Select the muscle relaxant size to be used in combination with the MFS anti-bruxism device taking the measurement from the mold of the patient's dental configuration.
2. Assemble the MFS anti-bruxism device, fitting the bite plate to the muscle relaxant's anterior groove.

3. Insert the medical device into the patient's mouth, placing the muscle relaxant at the buccal vestibule, outside the dental arches. The bite plate must be inserted between the lower and upper incisors.
4. A night-time regimen is indicated, although the patient may start wearing the device at home and during daytime, to get used to the device.

Instructions for Use: refer to page 0, fig. 7

4. RESULTS:

Clinical results are noticeable:

- a. Immediately in patients with bruxism
- b. After approximately six months in patients with deep anterior overbite and intense masticatory patterns.

5. SIZES

The MFS anti-bruxism device comes in six sizes, numbered from one to six. The smaller sizes are used in the pediatric population, while the larger sizes have been designed for use in adult patients.

There are no size specifications for this device, however, it must be underscored that the bite plate works passively, while the muscle relaxant works actively on cheek and lip muscles, pushing the muscles outwards (cheek) and forward (lips).

6. CONTRAINDICATIONS

Should not be used in dolichocephalic patients (long faced pattern) or in patients under five years of age.

7. WARNINGS


Bite plates should not be placed in direct contact with any heat source, otherwise the devices could be deformed and lose their physical characteristics.

8. OBSERVATIONS

It is convenient to separate both elements (muscle relaxant and bite plate) whenever the anti-bruxism device is not being used in order to maintain the relaxant's planar shape and ensure maximum effectiveness when it is inserted into the mouth again.

This device may be used in combination with fixed orthodontic appliances.

OPEN BITE

1. WHAT IS THE OPEN BITE TREATMENT DEVICE?

This medical device has been designed for the treatment of open bite malocclusion and performs the following functions:

- The lateral bite planes which the patient will bite on stimulate the masticatory muscles. If these muscles gain tone, the posterior teeth in occlusion will intrude and the vertical posterior dimension of occlusion will diminish.
- Control the forward thrust of the tongue with the device's anterior shield.

2. INDICATIONS

The open bite medical device is especially indicated in:

- Patients with anterior open bite, lax masticatory musculature and an advanced tongue that appears between the teeth.
- Patients with lax masticatory musculature (dolichocephalic) with a tendency towards a long faced pattern.

3. WHEN AND HOW TO USE THE OPEN BITE TREATMENT DEVICE

Guidelines to ensure that the MFS open bite treatment device is correctly used:

1. Select the open bite device size taking the measurement from the mold of the patient's dental configuration.
2. The device is placed within the mouth, behind the lower dental arch.
3. The lateral bite planes should cover inferior premolars and molars.

4. A night-time regimen is indicated. Initially, if the patient does not get used to wearing the device during the night, it may be used during daytime.
5. Clinical results will be noticeable after six to twelve months with closing of anterior bite malocclusion.
6. This device may be used in combination with fixed orthodontic appliances.

Instructions for Use: refer to page 0, fig. 8

4. SIZES

The open bite device comes in six sizes, numbered from one to six. The smaller sizes are used in the pediatric population, while the larger sizes have been designed for use in adult patients.

This medical device works passively and exerts no force, and this determines which size to use.

5. CONTRAINDICATIONS

Should not be used in brachycephalic patients (broad faced pattern).

6. WARNINGS

Should not be placed in direct contact with any heat source, otherwise the devices could be deformed and lose their physical characteristics.

7. OBSERVATIONS

In patients with limited mobility or with hypertrophic tonsils, we recommend addressing those conditions before using the device to ensure optimal performance.

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